


Anti-Taz Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Gzv5

Item Number: CM0020186



Introduction

High-quality rabbit polyclonal antibody against human TAZ (WWTR1). Validated for WB, IF/ICC, and ELISA; cross-reacts with human and mouse. Ideal for Hippo signaling and transcription regulation research. - Q9GZV5.

[Learn More](#)

Parameter	Value
Product Name	TAZ Rabbit pAb
Remarks/Alias	TAZ
Species	Human
Gene ID (Human)	25937
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 300-400 of human TAZ (NP_056287.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9GZV5
Protein Weight	44kDa
Shipping	Ice bag