


## Anti-Taz Rabbit Polymab® For Wb, Ihc-P, Ip, Elisa - Q9Gzv5

Item Number: CM0011677



### Introduction

Anti-TAZ rabbit recombinant monoclonal antibody (PolymAb®) for WB, IHC-P, IP, ELISA. Reacts with human, mouse, rat TAZ (WWTR1), a Hippo pathway transcriptional coactivator. Ideal for organ size control and tumor suppression research.

[Learn More](#)

| Parameter        | Value                                                                                                                               |
|------------------|-------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | TAZ Rabbit PolymAb®                                                                                                                 |
| Remarks/Alias    | TAZ                                                                                                                                 |
| Species          | Human                                                                                                                               |
| Gene ID (Human)  | 25937                                                                                                                               |
| Gene ID          | 25937                                                                                                                               |
| Immunogen        | Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human TAZ (NP_056287.1) |
| Source           | Rabbit                                                                                                                              |
| Category         | Monoclonal Antibodies                                                                                                               |
| Application      | WB, IHC-P, IP, ELISA                                                                                                                |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                   |
| SWISS            | Q9GZV5                                                                                                                              |
| Protein Weight   | 44kDa                                                                                                                               |
| Shipping         | Ice bag                                                                                                                             |