


Anti-Tbc1D1 Polyclonal Antibody For Wb And Elisa - Q86Ti0

Item Number: CM0025216



Introduction

TBC1D1 rabbit polyclonal antibody validated for Western blot and ELISA with human, mouse, and rat reactivity. Detects a 133 kDa band. Useful for investigating GLUT4 translocation, glucose uptake, and cell cycle regulation.

[Learn More](#)

Parameter	Value
Product Name	TBC1D1 Rabbit pAb
Remarks/Alias	TBC; TBC1; TBC1D1
Species	Human
GeneID (Human)	23216
GeneID	23216
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1069-1168 of human TBC1D1 (NP_055988.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86TI0
Protein Weight	133kDa
Shipping	Ice bag