


Anti-Tbc1D4 Polyclonal Antibody For Wb, If/Icc, Elisa - O60343

Item Number: CM0021901



Introduction

Anti-TBC1D4 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA. Recognizes human TBC1D4/AS160, cross-reactive with mouse and rat. Ideal for glucose metabolism and diabetes research.

[Learn More](#)

Parameter	Value
Product Name	TBC1D4 Rabbit pAb
Remarks/Alias	AS160; NIDDM5; TBC1D4
Species	Human
GeneID	9882
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1059-1298 of human TBC1D4 (NP_055647.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O60343
Protein Weight	147kDa
Shipping	Ice bag