

Anti-Tbcb Polyclonal Antibody For Wb, Ihc-P, Elisa - Q99426

Item Number: CM0015804

Introduction

Rabbit polyclonal antibody against human TBCB (tubulin-folding cofactor B). Validated for WB, IHC-P, ELISA, with cross-reactivity to mouse and rat. Targets a 27kDa protein involved in tubulin heterodimer assembly.

[Learn More](#)

Parameter	Value
Product Name	TBCB Rabbit pAb
Remarks/Alias	CG22; CKAP1; CKAPI; TBCB
Species	Human
GeneID (Human)	1155
GeneID	1155
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-244 of human TBCB (NP_001272.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99426
Protein Weight	27kDa
Shipping	Ice bag