

Anti-Tbcc Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q15814

Item Number: CM0023465

Introduction

Rabbit polyclonal antibody targeting TBCC (Tubulin-specific chaperone C), validated for WB, IF/ICC, ELISA. Recognizes human, mouse, rat. Suitable for studying tubulin folding and retinal biology.

[Learn More](#)

Parameter	Value
Product Name	TBCC Rabbit pAb
Remarks/Alias	CFC; TBCC
Species	Human
Gene ID (Human)	6903
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-346 of human TBCC (NP_003183.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15814
Protein Weight	39kDa
Shipping	Ice bag