

Anti-Tbccd1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Nvr7

Item Number: CM0019239

Introduction

Rabbit polyclonal antibody targeting human TBCCD1, validated for WB, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Suitable for investigating centrosome and Golgi positioning roles in cell shape and migration.

[Learn More](#)

Parameter	Value
Product Name	TBCCD1 Rabbit pAb
Remarks/Alias	TBCCD1
Species	Human
GeneID (Human)	55171
GeneID	55171
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 456-557 of human TBCCD1 (NP_060608.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NVR7
Protein Weight	64kDa
Shipping	Ice bag