

Anti-Tbce Polyclonal Antibody For If-P And Elisa - Q15813

Item Number: CM0022412



Introduction

Rabbit polyclonal antibody targeting human TBCE (tubulin-specific chaperone E), validated for immunofluorescence-paraffin and ELISA. Cross-reacts with mouse. Ideal for studying microtubule organization and neuronal maintenance.

[Learn More](#)

Parameter	Value
Product Name	TBCE Rabbit pAb
Remarks/Alias	HRD; KCS; KCS1; pac2; PEAMO; TBCE
Species	Human
GeneID (human)	6905
GeneID	6905
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 328-527 of human TBCE (NP_003184.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF-P, ELISA
Cross Reactivity	Mouse
SWISS	Q15813
Protein Weight	59kDa
Shipping	Ice bag