


Anti-Tbk1/Nak Monoclonal Antibody (Ko Validated) For Wb - Q9Uhd2

Item Number: CM0012383



Introduction

High-specificity rabbit anti-TBK1/NAK monoclonal antibody validated for WB, IF/ICC, IHC-P, IP, and ELISA. Reacts with human, mouse, and rat TBK1. Ideal for innate immunity, inflammation, and autophagy research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] TBK1/NAK Rabbit mAb
Remarks/Alias	NAK; T2K; IIAE8; FTDALS4; TBK1/NAK
Species	Human
GeneID (Human)	29110
GeneID	29110
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human TBK1/NAK (Q9UHD2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UHD2
Protein Weight	84 kDa
Shipping	Ice bag