


Anti-Tbk1/Nak Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Uhd2

Item Number: CM0030408



Introduction

High-quality rabbit polyclonal antibody against TBK1/NAK, validated for WB, IHC-P, ELISA with reactivity to human, mouse, rat. Ideal for innate immunity, inflammation, and autophagy research.

[Learn More](#)

Parameter	Specification
Product Name	TBK1/NAK Rabbit pAb
Remarks/Alias	NAK; T2K; IIAE8; FTDALS4; TBK1/NAK
Species	Human
GeneID	29110
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 500-600 of human TBK1/NAK (NP_037386.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UHD2
Protein Weight	84kDa
Shipping	Ice bag