

Anti-Tbl1X Polyclonal Antibody For Wb And Elisa - O60907

Item Number: CM0018274



Introduction

Rabbit polyclonal antibody against TBL1X (O60907) for WB and ELISA. Detects TBL1X in human, mouse, and rat samples. TBL1X is an F-box-like protein that recruits the proteasome to nuclear receptors, essential for transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	TBL1X Rabbit pAb
Remarks/Alias	EBI; TBL1; CHNG8; SMAP55; TBL1X
Species	Human
GeneID (Human)	6907
GeneID	6907
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60907
Protein Weight	62kDa
Shipping	Ice bag