

Anti-Tbl1Xr1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Bzk7

Item Number: CM0015664



Introduction

Rabbit polyclonal antibody against human TBL1XR1, validated for WB, IF/ICC, and ELISA. Reacts with human and mouse. Suitable for studying nuclear receptor-mediated transcription and ubiquitin-proteasome pathway.

[Learn More](#)

Parameter	Value
Product Name	TBL1XR1 Rabbit pAb
Remarks/Alias	C21; DC42; IRA1; MRD41; TBLR1; TBL1XR1
Species	Human
GeneID (Human)	79718
GeneID	79718
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-170 of human TBL1XR1 (NP_078941.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BZK7
Protein Weight	56kDa
Shipping	Ice bag