

Anti-Tbl2 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Y4P3

Item Number: CM0012974

Introduction

Anti-TBL2 Rabbit pAb validated for Western blot and ELISA, cross-reacts with human, mouse, and rat. Ideal for studying Williams-Beuren syndrome and WD repeat protein biology. Predicted molecular weight 50 kDa.

[Learn More](#)

Parameter	Value
Product Name	TBL2 Rabbit pAb
Remarks/Alias	WBSCR13; WS-betaTRP; TBL2
Species	Human
GeneID (Human)	26608
GeneID	26608
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-430 of human TBL2 (NP_036585.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y4P3
Protein Weight	50kDa
Shipping	Ice bag