

Anti-Tbp Rabbit Polyclonal Antibody For Wb And Elisa - P20226

Item Number: CM0023607



Introduction

High-quality anti-TBP rabbit polyclonal antibody validated for Western blot and ELISA applications. Detects human, mouse, and rat TATA-box-binding protein, a key transcription factor in RNA polymerase II initiation. Ideal for transcriptional regulation research.

[Learn More](#)

Parameter	Value
Product Name	TBP Rabbit pAb
Aliases	HDL4, TBP1, GTF2D, SCA17, TFIID, GTF2D1, TBP
Immunogen Species	Human
Gene ID (Human)	6908
Immunogen	Synthetic peptide within amino acids 1-100 of human TBP (NP_001165556.1)
Host Species	Rabbit
Antibody Type	Polyclonal Antibodies
Applications	WB, ELISA
Cross-Reactivity	Human, Mouse, Rat
UniProt Accession	P20226
Molecular Weight	38kDa
Shipping Method	Ice bag