

Anti-Tbr1 Monoclonal Antibody For Wb, If-P, Elisa - Q16650

Item Number: CM0003135

Introduction

Rabbit monoclonal antibody targeting human TBR1 (T-box brain protein 1), validated for WB, IF-P, and ELISA. Cross-reacts with mouse. Suitable for cortical development and neurodevelopmental disorder research.

[Learn More](#)

Parameter	Value
Product Name	TBR1 Rabbit mAb
Remarks/Alias	AUTS5; IDDAS; TBR-1; TES-56; TBR1
Species	Human
GeneID (human)	10716
GeneID	10716
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human TBR1 (Q16650)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse
SWISS	Q16650
Protein Weight	74 kDa
Shipping	Ice bag