

Anti-Tbrg4 Polyclonal Antibody For Wb - Q969Z0

Item Number: CM0023971



Introduction

CamelBio offers Anti-TBRG4 polyclonal antibody (Q969Z0) validated in WB, IF/ICC, ELISA; reacts with human, mouse, rat. Ideal for mitochondrial RNA processing and cell cycle research.

[Learn More](#)

Parameter	Details
Product Name	TBRG4 Rabbit pAb
Remarks/Alias	CPR2; FASTKD4; TBRG4
Species	Human
GeneID (human)	9238
GeneID	9238
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 382-631 of human TBRG4 (NP_004740.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q969Z0
Protein Weight	71kDa
Shipping	Ice bag