

Anti-Tbx19 Rabbit Polyclonal Antibody For Wb, Elisa - O60806

Item Number: CM0017819

Introduction

Rabbit polyclonal antibody against human TBX19 (TPIT), a nuclear T-box transcription factor. Validated for Western blot (WB) and ELISA. Ideal for studying pituitary development and transcriptional regulation.

[Learn More](#)

Parameter	Value
Product Name	TBX19 Rabbit pAb
Remarks/Alias	TPIT; TBS19; dj747L4.1; TBX19
Species	Human
GeneID	9095 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 229-448 of human TBX19 (NP_005140.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O60806
Protein Weight	48kDa
Shipping	Ice bag