

Anti-Tbx20 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Umr3

Item Number: CM0023327

Introduction

Rabbit polyclonal antibody against human TBX20, validated for WB, IHC-P, ELISA. Cross-reacts with human and rat. Suitable for cardiac development research and transcription factor studies.

[Learn More](#)

Parameter	Value
Product Name	TBX20 Rabbit pAb
Remarks/Alias	ASD4; TBX20
Species	Human
GeneID (Human)	57057
GeneID	57057
Immunogen	Recombinant Protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9UMR3
Protein Weight	49kDa
Shipping	Ice bag