

Anti-Tbx21 Rabbit Polyclonal Antibody For Wb, Elisa - Q9UI17

Item Number: CM0028528

Introduction

TBX21 Rabbit pAb is a polyclonal antibody raised against the C-terminal region of human T-bet, validated for WB and ELISA. Recognizes human, mouse, and rat TBX21, a master transcription factor for Th1 cell development, crucial for immunology and autoimmune disease research.

[Learn More](#)

Parameter	Value
Product Name	TBX21 Rabbit pAb
Remarks/Alias	TBET; IMD88; T-PET; T-bet; TBLYM; TBX21
Species	Human
Gene ID (Human)	30009
Gene ID	30009
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UL17
Protein Weight	58kDa
Shipping	Ice bag