

Anti-Tbx22 Polyclonal Antibody For Wb, Elisa - Q9Y458

Item Number: CM0026114

Introduction

Rabbit polyclonal antibody targeting human TBX22 (Q9Y458). Suitable for WB and ELISA, with cross-reactivity to human, mouse, and rat. Ideal for studying palatogenesis and transcriptional regulation in developmental biology.

[Learn More](#)

Parameter	Value
Product Name	TBX22 Rabbit pAb
Remarks/Alias	CPX; CLPA; TBXX; ABERS; dj795G23.1; TBX22
Species	Human
GeneID (Human)	50945
GeneID	50945
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 291-520 of human TBX22 (NP_058650.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y458
Protein Weight	58kDa
Shipping	Ice bag