

Anti-Tbx6 Rabbit Polyclonal Antibody For Wb, Elisa - 095947

Item Number: CM0020574

Introduction

Affinity-purified rabbit polyclonal antibody against human TBX6 (T-box transcription factor), validated for Western blot (WB) and ELISA, targeting mesoderm and left-right axis research applications.

[Learn More](#)

Parameter	Value
Product Name	TBX6 Rabbit pAb
Remarks/Alias	SCDO5; TBX6
Species	Human
GeneID (Human)	6911
GeneID	6911
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human TBX6 (NP_004599.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O95947
Protein Weight	47kDa
Shipping	Ice bag