


Anti-Tc10/Rhoq Monoclonal Antibody For Wb, Ihc-P, Elisa - P17081

Item Number: CM0002752



Introduction

Rabbit monoclonal antibody targeting TC10/RHOQ (RhoQ) small GTPase. Validated for WB, IHC-P, ELISA with human, mouse, rat cross-reactivity. Ideal for studying cell polarization, filopodia formation, and CFTR trafficking.

[Learn More](#)

Parameter	Value
Product Name	TC10/RHOQ Rabbit mAb
Remarks/Alias	ARHQ; TC10; TC10A; RASL7A; HEL-S-42; TC10/RHOQ
Species	Human
GeneID	23433
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-205 of human TC10/RHOQ (P17081).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17081
Protein Weight	23kDa
Shipping	Ice bag