

Anti-Tcap Polyclonal Antibody For Wb, Elisa - 015273

Item Number: CM0014665



Introduction

Rabbit polyclonal antibody against human TCAP (telethonin). Validated for Western blot and ELISA, with cross-reactivity to mouse and rat. Suitable for muscle sarcomere assembly and disease research.

[Learn More](#)

Parameter	Value
Product Name	TCAP Rabbit pAb
Remarks/Alias	TELE; CMD1N; CMH25; T-cap; LGMD2G; LGMDR7; telethonin; TCAP
Species	Human
Gene ID (Human)	8557
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-167 of human TCAP (NP_003664.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	015273
Protein Weight	19kDa
Shipping	Ice bag