

Anti-Tcap Rabbit Monoclonal Antibody For Wb, If-P, Elisa - O15273

Item Number: CM0004897



Introduction

High-quality Anti-TCAP rabbit monoclonal antibody for WB, IF-P, ELISA. Detects human, mouse, rat telethonin. Ideal for sarcomere research and muscle disease studies.

[Learn More](#)

Parameter	Value
Product Name	TCAP Rabbit mAb
Remarks/Alias	TELE; CMD1N; CMH25; T-cap; LGMD2G; LGMDR7; telethonin; TCAP
Species	Human
Gene ID (human)	8557
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-167 of human TCAP (NP_003664.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15273
Protein Weight	19kDa
Shipping	Ice bag