

Anti-Tcap Rabbit Polyclonal Antibody For Wb, Elisa - O15273

Item Number: CM0025937



Introduction

High-quality rabbit polyclonal antibody against human telethonin (TCAP). Validated for Western blotting (WB) and ELISA. Suitable for studying sarcomere assembly, cardiomyopathy, and limb-girdle muscular dystrophy. Reacts with human samples.

[Learn More](#)

Parameter	Value
Product Name	TCAP Rabbit pAb
Remarks/Alias	TELE; CMD1N; CMH25; T-cap; LGMD2G; LGMDR7; telethonin; TCAP
Species	Human
Gene ID (Human)	8557
Gene ID	8557
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-167 of human TCAP (NP_003664.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O15273
Protein Weight	19kDa
Shipping	Ice bag