

Anti-Tcea1 Monoclonal Antibody For Wb, If/Icc, Elisa - P23193

Item Number: CM0002279

Introduction

Rabbit monoclonal antibody targeting human TCEA1, a transcription elongation factor. Validated for WB, IF/ICC, ELISA with cross-reactivity to mouse and rat. Ideal for transcription and gene regulation research.

[Learn More](#)

Parameter	Value
Product Name	TCEA1 Rabbit mAb
Remarks/Alias	SII; TCEA; TF2S; GTF2S; TFIIS; TCEA1
Species	Human
GeneID (Human)	6917
GeneID	6917
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human TCEA1 (P23193)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23193
Protein Weight	34kDa
Shipping	Ice bag