


Anti-Tceb1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q15369

Item Number: CM0017069



Introduction

Rabbit polyclonal antibody targeting human TCEB1 (Elongin-C). Validated for WB, IHC-P, ELISA, and cross-reacts with mouse and rat. Ideal for studies on transcription elongation and ubiquitin ligase complexes.

[Learn More](#)

Parameter	Value
Product Name	TCEB1 Rabbit pAb
Remarks/Alias	SIII; TCEB1
Species	Human
GeneID (Human)	6921
GeneID	6921
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-112 of human TCEB1 (NP_001191791.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15369
Protein Weight	12kDa
Shipping	Ice bag