

Anti-Tceb2 Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q15370

Item Number: CM0015197



Introduction

CamelBio offers Anti-TCEB2 Polyclonal Antibody validated for WB, IHC-P, IP, ELISA. Targets human Elongin-B, a transcription elongation factor and ubiquitin ligase component. Ideal for gene regulation and cancer research.

[Learn More](#)

Parameter	Specification
Product Name	TCEB2 Rabbit pAb
Remarks/Alias	SIII; TCEB2
Species	Human
Gene ID	6923
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-118 of human TCEB2 (NP_009039.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human
SWISS	Q15370
Protein Weight	13kDa
Shipping	Ice bag