

Anti-Tceb2 Rabbit Monoclonal Antibody For Wb And Elisa - Q15370

Item Number: CM0006760



Introduction

CamelBio supplies anti-TCEB2 rabbit monoclonal antibody (clone unspecified) for WB and ELISA, with cross-reactivity to human, mouse, and rat. Essential for transcription elongation and ubiquitination research in IVD development.

[Learn More](#)

Parameter	Value
Product Name	TCEB2 Rabbit mAb
Remarks/Alias	SIII; TCEB2
Species	Human
GeneID (Human)	6923
GeneID	6923
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15370
Protein Weight	13kDa
Shipping	Ice bag