


Anti-Tcerg1 Polyclonal Antibody For Wb, Elisa - O14776

Item Number: CM0014995



Introduction

Rabbit polyclonal antibody targeting human TCERG1 (CA150/TAF2S), validated for Western blot and ELISA. Ideal for studying transcription elongation regulation and HIV-1 Tat-mediated activation.

[Learn More](#)

Parameter	Value
Product Name	TCERG1 Rabbit pAb
Remarks/Alias	Urn1; CA150; TAF2S; TCERG1
Species	Human
Gene ID (Human)	10915
Gene ID	10915
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 550-650 of human TCERG1 (NP_006697.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O14776
Protein Weight	124kDa
Shipping	Ice bag