

Anti-Tcf12 Rabbit Polyclonal Antibody For Wb, Elisa - Q99081

Item Number: CM0026748



Introduction

Rabbit polyclonal antibody against human TCF12 (Transcription factor 12), suitable for WB and ELISA. Recognizes human, mouse, and rat TCF12. Ideal for studying neuronal differentiation and transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	TCF12 Rabbit pAb
Remarks/Alias	HEB; p64; CRS3; HH26; HTF4; TCF-12; bHLHb20; HsT17266; TCF12
Species	Human
Gene ID (Human)	6938
Gene ID	6938
Immunogen	Recombinant fusion protein containing amino acids 443-682 of human TCF12
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99081
Protein Weight	73kDa
Shipping	Ice bag