

Anti-Tcf19 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y242

Item Number: CM0030988

Introduction

Rabbit polyclonal antibody against human TCF19 (transcription factor 19). Validated for WB and ELISA. Cross-reacts with human, mouse, and rat. Ideal for cell cycle and transcription regulation research.

[Learn More](#)

Parameter	Value
Product Name	TCF19 Rabbit pAb
Remarks/Alias	SC1; TCF-19; TCF19
Species	Human
GeneID (Human)	6941
GeneID	6941
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-290 of human TCF19 (NP_009040.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y242
Protein Weight	37kDa
Shipping	Ice bag