

Anti-Tcf1/Tcf7 Monoclonal Antibody For Wb, Elisa - P36402

Item Number: CM0002374

Introduction

Rabbit monoclonal antibody against human TCF1/TCF7 (TCF-1, P36402). Validated for WB, ELISA; cross-reacts with mouse and rat. Ideal for studying T-cell differentiation and Wnt signaling.

[Learn More](#)

Parameter	Description
Product Name	TCF1/TCF7 Rabbit mAb
Remarks/Alias	TCF-1; TCF1/TCF7
Species	Human
GeneID (Human)	6932
GeneID	6932
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human TCF1/TCF7 (NP_003193.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P36402
Protein Weight	42kDa
Shipping	Ice bag