

Anti-Tcf1/Tcf7 Rabbit Polyclonal Antibody For Wb, Elisa - P36402

Item Number: CM0029034

Introduction

High-quality rabbit polyclonal antibody against TCF1/TCF7, validated for WB and ELISA. Detects human, mouse, and rat TCF7, a transcription factor critical for T-cell development. Ideal for immunodetection and research applications.

[Learn More](#)

Parameter	Value
Product Name	TCF1/TCF7 Rabbit pAb
Remarks/Alias	TCF-1; TCF1/TCF7
Species	Human
GenelD (human)	6932
GenelD	6932
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human TCF1/TCF7 (NP_003193.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P36402
Protein Weight	42kDa
Shipping	Ice bag