

Anti-Tcf21 (Pod1/Bhlha23) Rabbit Polyclonal Antibody For Wb And Elisa - O43680

Item Number: CM0018823

Introduction

Rabbit polyclonal antibody against TCF21 (POD1/bHLHa23), validated for WB and ELISA. Detects human and mouse TCF21, a transcription factor involved in kidney and lung development.

[Learn More](#)

Parameter	Value
Product Name	TCF21 Rabbit pAb
Remarks/Alias	POD1; bHLHa23; TCF21
Species	Human
Gene ID (Human)	6943
Gene ID	6943
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-100 of human TCF21 (NP_003197.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O43680
Protein Weight	20kDa
Shipping	Ice bag