


Anti-Tcf3 Rabbit Polyclonal Antibody For Wb, Elisa - P15923

Item Number: CM0023946



Introduction

TCF3 Rabbit pAb is a polyclonal antibody validated for WB and ELISA, recognizing human, mouse, and rat TCF3 (68 kDa), a transcription factor involved in B-cell differentiation and neuronal development.

[Learn More](#)

Parameter	Value
Product Name	TCF3 Rabbit pAb
Remarks/Alias	E2A; E47; p75; AGM8; ITF1; VDIR; AGM8A; AGM8B; TCF-3; bHLHb21; TCF3
Species	Human
GeneID (Human)	6929
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 450-525 of human TCF3 (NP_003191.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15923
Protein Weight	68kDa
Shipping	Ice bag