


Anti-Tcf4 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P15884

Item Number: CM0022658



Introduction

High-quality anti-TCF4 rabbit polyclonal antibody for applications including WB, IHC-P, IF/ICC, and ELISA. Ideal for studying transcriptional regulation, neuronal differentiation, and TCF4-associated diseases.

[Learn More](#)

Parameter	Value
Product Name	TCF4 Rabbit pAb
Remarks/Alias	E2-2; FCD2; ITF2; PTHS; SEF2; CDG2T; FECD3; ITF-2; SEF-2; TCF-4; SEF2-1; SEF2-1A; SEF2-1B; SEF2-1D; bHLHb19; TCF4
Species	Human
Gene ID (Human)	6925
Gene ID	6925
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human TCF4 (NP_003190.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15884
Protein Weight	71kDa
Shipping	Ice bag