


Anti-Tcf4/Tcf7L2 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Nqb0

Item Number: CM0002678



Introduction

Recombinant rabbit monoclonal antibody targeting human TCF4/TCF7L2, a Wnt signaling transcription factor. Validated for Western blot, IHC-P, and ELISA. Supports cancer and stem cell research.

[Learn More](#)

Parameter	Value
Product Name	TCF4/TCF7L2 Rabbit mAb
Remarks/Alias	TCF4; TCF-4; TCF4/TCF7L2
Species	Human
GeneID (human)	6934
GeneID	6934
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human TCF4/TCF7L2 (Q9NQB0)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q9NQB0
Protein Weight	68 kDa
Shipping	Ice bag