

Anti-Tcf4/Tcf7L2 Polyclonal Antibody For Wb, Elisa - Q9Nqb0

Item Number: CM0022661



Introduction

Rabbit polyclonal antibody targeting TCF4/TCF7L2 (Transcription factor 7-like 2) for Western blot and ELISA applications. Recognizes human and mouse TCF4/TCF7L2, a Wnt signaling transcription factor with roles in cancer and epithelial stem cell maintenance.

[Learn More](#)

Parameter	Value
Product Name	TCF4/TCF7L2 Rabbit pAb
Remarks/Alias	TCF4; TCF-4; TCF4/TCF7L2
Species	Human
GeneID (Human)	6934
GeneID	6934
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 124-188 of human TCF4/TCF7L2 (NP_001354872.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NQB0
Protein Weight	68kDa
Shipping	Ice bag