

Anti-Tcf7L1 Polyclonal Antibody For Wb And Elisa - Q9Hcs4

Item Number: CM0026326

Introduction

Rabbit polyclonal antibody against TCF7L1 (TCF-3) validated for Western blot and ELISA. Reacts with human, mouse, and rat. Useful for Wnt signaling and epidermal differentiation studies.

[Learn More](#)

Parameter	Value
Product Name	TCF7L1 Rabbit pAb
Remarks/Alias	TCF3; TCF-3; TCF7L1
Species	Human
GeneID (Human)	83439
GeneID	83439
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 260-340 of human TCF7L1 (NP_112573.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human,Mouse,Rat
SWISS	Q9HCS4
Protein Weight	63kDa
Shipping	Ice bag