

Anti-Tcirg1 Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q13488

Item Number: CM0017664

Introduction

Rabbit polyclonal antibody against TCIRG1 (V-ATPase subunit a3), validated for WB, IF/ICC, IHC-P, ELISA. Reacts with human, mouse, rat. Suitable for osteoclast and T-cell research.

[Learn More](#)

Parameter	Value
Product Name	TCIRG1 Rabbit pAb
Remarks/Alias	a3; Stv1; Vph1; Atp6i; OC116; OPTB1; TIRC7; ATP6N1C; ATP6V0A3; OC-116kDa; TCIRG1
Species	Human
GeneID (Human)	10312
GeneID	10312
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human TCIRG1 (NP_006010.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13488
Protein Weight	93 kDa
Shipping	Ice bag