

Anti-Tcirg1 Rabbit Monoclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - Q13488

Item Number: CM0006963

Introduction

Rabbit monoclonal antibody targeting human TCIRG1 (V-ATPase subunit a3). Validated for WB, IF, IHC, ELISA. Reacts with human, mouse, rat. Useful for osteoclast, immune regulation, and pH homeostasis studies.

[Learn More](#)

Parameter	Value
Product Name	TCIRG1 Rabbit mAb
Remarks/Alias	a3; Stv1; Vph1; Atp6i; OC116; OPTB1; TIRC7; ATP6N1C; ATP6V0A3; OC-116kDa; TCIRG1
Species	Human
Gene ID (Human)	10312
Gene ID	10312
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-356 of human TCIRG1 (NP_006010.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13488
Protein Weight	93kDa
Shipping	Ice bag