

Anti-Tcl1 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P56279

Item Number: CM0008045



Introduction

Anti-TCL1 rabbit monoclonal antibody validated for WB, IF/ICC, ELISA. Targets human TCL1 (13 kDa), enhancing AKT phosphorylation and promoting cell survival. Suitable for leukemia/lymphoma research.

[Learn More](#)

Parameter	Value
Product Name	TCL1 Rabbit mAb
Remarks/Alias	TCL1
Species	Human
Gene ID (Human)	8115
Gene ID	8115
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-114 of human TCL1 (NP_068801.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P56279
Protein Weight	13kDa
Shipping	Ice bag