

Anti-Tcl1A Polyclonal Antibody For Wb, Elisa - P56279

Item Number: CM0017057



Introduction

TCL1A Rabbit pAb validated for WB and ELISA in human samples. Targets the T-cell leukemia/lymphoma protein 1A, an AKT co-activator involved in cell proliferation and survival. Ideal for oncology and lymphoid research.

[Learn More](#)

Parameter	Value
Product Name	TCL1A Rabbit pAb
Remarks/Alias	TCL1; TCL1A
Species	Human
GeneID (Human)	8115
GeneID	8115
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-114 of human TCL1A (NP_001092195.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P56279
Protein Weight	13kDa
Shipping	Ice bag