

Anti-Tcl1B Rabbit Polyclonal Antibody For Wb, Elisa - 095988

Item Number: CM0016925

Introduction

High-quality anti-TCL1B rabbit polyclonal antibody validated for Western blotting and ELISA. Detects human, mouse, rat TCL1B (predicted ~15 kDa). Useful for studying proto-oncogene signaling and AKT activation pathways.

[Learn More](#)

Parameter	Value
Product Name	TCL1B Rabbit pAb
Remarks/Alias	TML1; SYN-1; TCL1B
Species	Human
Gene ID	9623
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-128 of human TCL1B (NP_004909.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	095988
Protein Weight	15kDa
Shipping	Ice bag