

Anti-Tcn2 (Transcobalamin-2) Rabbit Polyclonal Antibody For Wb, Elisa - P20062

Item Number: CM0029996



Introduction

Rabbit polyclonal antibody targeting human Transcobalamin-2 (TCN2), suitable for Western blotting and ELISA, with cross-reactivity to mouse TCN2. Suitable for cobalamin transport and metabolism research.

[Learn More](#)

Parameter	Value
Product Name	TCN2 Rabbit pAb
Remarks/Alias	II; TC; TC2; TC-2; TCII; TC II; D22S676; D22S750; TCN2
Species	Human
Gene ID	6948 (Human)
Immunogen	Recombinant fusion protein containing amino acids 19-250 of human TCN2 (NP_001171655.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P20062
Protein Weight	48kDa
Shipping	Ice bag