

Anti-Tcof1 Polyclonal Antibody For Wb And Elisa - Q13428

Item Number: CM0030197

Introduction

Rabbit polyclonal antibody against human TCOF1 (Treacle protein), suitable for Western blot and ELISA. Supports research in nucleolar function, RNA polymerase I regulation, and neural crest specification.

[Learn More](#)

Parameter	Value
Product Name	TCOF1 Rabbit pAb
Remarks/Alias	TCS; MFD1; TCS1; treacle; TCOF1
Species	Human
GeneID (Human)	6949
GeneID	6949
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 659-958 of human TCOF1 (NP_001008657.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q13428
Protein Weight	152kDa
Shipping	Ice bag