


Anti-Tcp1 Alpha Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P17987

Item Number: CM0007751



Introduction

CamelBio's anti-TCP1 alpha rabbit monoclonal antibody (P17987) validated for WB, IHC-P, ELISA. Detects human TCP1 alpha, a TRiC chaperonin subunit involved in protein folding and cytoskeleton assembly.

[Learn More](#)

Parameter	Value
Product Name	TCP1 alpha Rabbit mAb
Remarks/Alias	CCT1; CCTa; D6S230E; CCT-alpha; TCP-1-alpha; TCP1 alpha
Species	Human
GeneID	6950
Immunogen	Synthetic Peptide (A synthetic peptide corresponding to a sequence within amino acids 1-100 of human TCP1 alpha (P17987))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P17987
Protein Weight	60kDa
Shipping	Ice bag