


Anti-Tcp1 Beta/Cct2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P78371

Item Number: CM0005703



Introduction

Rabbit monoclonal antibody targeting human TCP1 beta/CCT2 (P78371), a TRiC chaperonin subunit. Validated for WB, IHC-P, ELISA in human, mouse, rat samples. Ideal for folding and ciliogenesis research.

[Learn More](#)

Parameter	Value
Product Name	TCP1 beta/CCT2 Rabbit mAb
Remarks/Alias	CCTB; 99D8.1; PRO1633; CCT-beta; HEL-S-100n; TCP-1-beta; TCP1 beta/CCT2
Species	Human
GeneID (Human)	10576
GeneID	10576
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 436-535 of human TCP1 beta/CCT2 (P78371).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P78371
Protein Weight	57kDa
Shipping	Ice bag