


Anti-Tcp1 Beta/Cct2 Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P78371

Item Number: CM0028257



Introduction

Rabbit polyclonal antibody targeting human TCP1 beta/CCT2 (CCT2), validated for Western blot, immunofluorescence, immunoprecipitation, and ELISA. Detects both human and mouse TCP1 beta. Ideal for chaperonin complex and protein folding research.

[Learn More](#)

Parameter	Value
Product Name	TCP1 beta/CCT2 Rabbit pAb
Remarks/Alias	CCTB; 99D8.1; PRO1633; CCT-beta; HEL-S-100n; TCP-1-beta; TCP1 beta/CCT2
Species	Human
GeneID	10576
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 206-535 of human TCP1 beta/TCP1 beta/CCT2 (NP_006422.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P78371
Protein Weight	57kDa
Shipping	Ice bag